

prietest™ R-SYS...Clinical Chemistry Reagents

TOTAL CHOLESTEROL TEST KIT

INTENDED USE :

- Quantitative in vitro determination of concentration of Total cholesterol in serum or plasma on photometric systems.
- In vitro diagnostic test kit, for laboratory and professional use.
- This manual contains instructions for operation by qualified personnel only.

ORDERING INFORMATION

Pack Size
6 X 26 ml

Cat No.
CHO0626RS

CLINICAL SIGNIFICANCE : Cholesterol is both coming from food and synthesized by the human body, mainly in hepatic and intestinal cells. Cholesterol is a component of cell membranes and a precursor for steroid hormones and bile acids synthesized by body cells and absorbed with food. Cholesterol is transported in plasma via lipoproteins, namely complexes between lipids and apolipoproteins. There are four classes of lipoproteins: high density lipoproteins (HDL), low density lipoproteins (LDL), very low density lipoproteins (VLDL) and chylomicrons. The determination of the individual total cholesterol (TC) level is used for screening purposes while for a better risk assessment it is necessary to measure additionally HDL-Cholesterol and LDL-Cholesterol.

METHOD : "CHOD-PAP": enzymatic photometric test, Trinder, End Point.

PRINCIPLE : Determination of cholesterol after enzymatic hydrolysis and oxidation. The colorimetric indicator is quinonimine which is generated from 4-aminoantipyrine and phenol by hydrogen peroxide under the catalytic action of peroxidase (Trinder's reaction).



REAGENTS :

COMPONENTS AND CONCENTRATIONS:

Pipes buffer	: 100 mmol/l
Cholesterol Oxidase	: > 100 U/l
Peroxidase	: > 500 U/l
Cholesterol Esterase	: > 150 U/l
4 - Amino Antipyrine	: 0.5 mmol/l
Phenol	: 10 mmol/l
Preservative & Stabilizer	

STORAGE INSTRUCTIONS AND REAGENT STABILITY : The reagent is stable up to the end of the indicated date of expiry on the vial label, if stored at 2 to 8°C, protected from light and contamination is avoided. Do not freeze the reagents.

The standard is stable up to the end of the indicated date of expiry on the vial label, if stored at 2 to 8°C.

Note: It has to be mentioned, that the measurement is not influenced by occasionally occurring color changes, as long as the absorbance of the reagent is < 0.3 at 510 nm.

WARNINGS AND PRECAUTIONS :

- The reagent contains sodium azide (0.95 g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
- Take the necessary precautions for the use of laboratory reagents.

WASTE MANAGEMENT : Please refer to local regulation requirements.

REAGENT PREPARATION : The reagent is ready to use.

MATERIALS REQUIRED BUT NOT PROVIDED : NaCl solution 9 g/l, General laboratory equipment, Analyser / Photometer, Pipettes etc.

SPECIMEN : Serum, heparinized plasma or EDTA plasma from fasting patients.

Stability : 7 days at 4 to 8°C, 3 months at -20°C

Discard contaminated specimens.

CONVERSION FACTOR : Cholesterol [mg/dl] X 0.0259 = Cholesterol [mmol/l]

CALIBRATION : For the calibration of automated photometric systems use of the commercially available calibrator is recommended.

QUALITY CONTROL : To ensure adequate quality, use of the commercially available priet serum is recommended.

PERFORMANCE CHARACTERISTICS :

MEASURING RANGE : The test has been developed to determine Cholesterol concentrations within a measuring range from 4 to 700 mg/dl. (0.104 to 18.13 mmol/L). When values exceed higher limit of the range, such samples should be diluted 1 + 1 with NaCl solution (9 g/l) and the result multiplied by 2.

SPECIFICITY / INTERFERENCES : No interference was observed by ascorbic acid up to 5 mg/dl, (283.9 µmol/L), bilirubin up to 20 mg/dl (342.00 µmol/L), hemoglobin up to 0.2 g/dl (2 g/L) and triglycerides up to 600 mg/dl (6.84 mmol/L). A list of drugs and other interfering substances with Cholesterol determination has been reported by Young et al.

SENSITIVITY / LIMIT OF DETECTION : The lower limit of detection is 4 mg/dl (0.104 mmol/L).

PRECISION :

Intra-assay precision n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
Sample 1	234.26	2.29	0.98
Sample 2	196.8	1.67	0.85
Sample 3	188.8	1.31	0.69
Inter-assay precision n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
Sample 1	145.8	2.43	1.66
Sample 2	91.89	1.71	1.86
Sample 3	259.9	4.66	1.79

METHOD COMPARISON : A comparison between Robonik Prietest Cholesterol (y) and a commercially available test (x) using 20 samples gave following results :

Linear Regression : y = 0.9602 x + 7.3582 mg/dl

Correlation Coefficient : r = 0.9762

REFERENCE RANGE :

Desirable	: ≤ 200 mg/dl	(≤ 5.18 mmol/L)
Borderline high risk	: 200 to 240 mg/dl	(5.18 to 6.22 mmol/L)
High risk	: > 240 mg/dl	(> 6.22 mmol/L)

It is recommended that each laboratory should assign its own normal range.

LITERATURE :

- Naito, H. K., *Coronary artery disease and disorders of lipid metabolism*, Clinical Chemistry: Theory, Analysis, Correlation, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S. C., (Mosby Inc. eds St Louis USA), (2003).
- Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed., (W.B. Saunders eds. Philadelphia USA),.
- Young DS. *Effects of drugs on Clinical Lab. Tests*, 4th ed. AACC Press, 1995.

INSTRUMENT APPLICATION : *autora / smara*

Reagent

CHOLESTEROL_ROBONIK

ALBUMIN_ROBONIK
ALKALINE PHOSPHATASE
ALKALINE PHOSPHATASE_R
AMYLASE_ROBONIK
ASO
ASO_ROBONIK
CALCIUM OCPC_ROBONIK
CALCIUM_ROBONIK
CHLORIDE_ROBONIK
CHOLESTEROL_ROBONIK
CREATININE_ROBONIK
CRP
CRP_ROBONIK
DIRECT BILIRUBIN_ROBONIK

Reagent Name :
Company :

CHOLESTEROL
ROBONIK

Mode Details :

Mode :
Filter :
Primary
Secondary

End Point
510
None

Volume Info

Sample Volume :
R1 Test Volume :
R2 Test Volume :

3
300
0

µl
µl
µl

Initial Absorbance :

Initial Abs.

Procedure

R1R2 in Same Cuvette
R1R2 in Different Cuvette

R1 + R2 + Sample
R1 + Sample + R2

R1 + Sample + R2

Reagent Blank
Sample Blank

Final Reading

R2-R1
R2/R1

New

Save

Delete

Cancel

Close

Reaction Time (sec)

Lag Time (T1) :
Lag Time (T2) :
Read Time (T3) :

600
0
0

s
s
s

Unit :

K Factor :

Blank
Factor
Single Calibrator
Multi Test Calibrator *

mg/dl
1

Single Calibrator
Multi Test Calibrator *

Calibrator Info :
No. Of Calibrators

1

*** Refer Pack Insert of Multi Test Calibrator to Input Value**

prietest is the Trade Mark of ROBONIK (INDIA) PVT.LTD., for *Clinical Chemistry Reagents*.



An ISO 9001 : 2015 Certified Company
An ISO 13485 : 2016 Certified Company

Manufactured and Marketed by:

ROBONIK (INDIA) PVT. LTD.,
Plot No. 3 & 4, MIDC Industrial Area, Morivali, Near Ladi Naka,
Ambarnath (West) - 421 501, District Thane, Maharashtra, INDIA,
Tel.: +91 (251) 2689000,
Email : customercare@robonikindia.com, Website : www.robonik.in

Toll Free No. 1800 5727 977



For in vitro diagnostic use	Catalogue Number
Consult Instructions for use	Lot No.
Manufacturer's Address	Date of Manufacture
Store at	Exp. Date